Manufacturer: An Hui Su Bo Lun Clothing Co., Ltd. Model Tested: SU KN95 Date Tested: April 22, 2020

These findings pertain to the An Hui Su Bo Lun Clothing Co., Ltd., SU KN95. The labeling for this product indicates that it meets GB2626-2006 KN95, the Chinese standard for Respiratory Protective Devices – Filtering Half Masks to Protect Against Particles – Requirements, Testing, Marking.

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found <u>here</u>.

A certificate of approval was not provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 95.10% and 28.45%, respectively. Nine of the respirators measured less than 95%. Only one respirator measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies (during known</u> <u>shortages)</u>.

## **Evaluation of International Respirators**

Test: Modified TEB-APR-STP-0059

Date Tested: April 22, 2020

Report Prepared: April 22, 2020

Manufacturer: An Hui Su Bo Lun Clothing Co., Ltd.

Item Tested: SU KN95

#### Country of Certification: China (GB2626-2006)

**Initial Filter Initial Percent** Maximum Flow Rate Filter Resistance Leakage Percent Leakage Filter Efficiency (Lpm) (mmH<sub>2</sub>O) (%) (%) 85 57.90 57.90 42.10 1 6.9 2 85 6.0 63.90 63.90 36.10 3 85 10.9 5.27 5.27 94.73 4 85 6.2 68.30 68.30 31.70 5 85 11.0 9.64 9.64 90.36 6 85 5.1 71.55 71.55 28.45 7 85 8.2 51.60 51.60 48.40 8 85 10.8 9.34 9.34 90.66 9 85 7.8 54.35 54.35 45.65 10 85 10.8 4.90 4.90 95.10 Minimum Filter Efficiency: 28.45 Maximum Filter Efficiency: 95.10

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NATIONAL Personal Protective Technology Laboratory

Pictures have been added to the end of this report.



陵用设服、 佩戴说明 面向口

1	偏影说明,	使用说明:
	[ ILay]	进入污染区之前请佩戴好口罩并检查其
	1. 面向口罩无鼻头	
	使鼻夹位于口罩 2. 将口罩抵住下巴	
	头带拉过头顶,	署千硚后 即更换
	耳朵下放。 3. 将上方头带拉边	本口罩不得改装、清洗 生头顶,置 本口罩不产生氧气,不能在氧气浓度低
	5. 将上方关带拉拉 于颈后耳朵上方	
	4. 将双手中指置于	金属鼻夹 区域中使用
	▲ 「 Y→Y」 中部、从中向保 梁形状向脸部招	
	来形状间腔部的 将其完全按压成	自动型山 仁志事员: 本口阜个能防护有害气体相差
	为止	海来形状 汽;本口罩不产生氧气,不能用于缺氧环境;不可用于存在油雾
	5. 在进入工作区域	1 战利尔克; 个时用于仔仕油萎的环境。
	须检查口罩与脸	部的密合
	检查方法: 1)用双手招	2注口面
	大力呼吸、如空	
	新调整鼻夹;如	
F 4	罩边缘溢出,应 头带,如不能取	
. 6	密合,应重复步	聚1-4 山山,中等小师
24	2) 如没有感	一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一一
7	出,可进入工作	区工作 电话: 400-0556-089
Conception of the local division of the loca	a constructed of the state of the	

产品合格证 产品名称: 折叠式口罩 执行标准: GB2626-2006 数 量: 1只、2只、5只、10只 生产日期: 生产日期: 2020年 03月 19日 检验日期: 2020年 03月 19日 检验员:01号 安徽苏博伦服饰有限公司 KN95 0B2626-2006



